### Next Generation Large-Format Spectroscopic Arrays for SOFIA

How a new generation of heterodyne arrays will justify continuing SOFIA's operation

Paul F. Goldsmith
Jet Propulsion Laboratory

**SOFIA Instrument Roadmap workshop II** 

July 27, 2020

With special thanks to Jose Siles, Jon Kawamura, and Imran Mehdi




## The Importance of Velocity-Resolved Spectroscopic Imaging


- Spectroscopy has for years been recognized as SOFIA's strongest area of contribution to astrophysics
  - Access to critical lines blocked by the atmosphere (fine structure lines including [CII], [NII], [OI]), unique molecular lines (HF), and high-J CO lines
- High resolution spectroscopy is even more favorable for SOFIA
  - Residual warm atmosphere and warm telescope DO NOT significantly degrade sensitivity of observations (this is not the case for photometry and low-resolution spectroscopy)
- The value of such data for unraveling key processes related to star formation, the life cycle of the ISM, and feedback is enormously enhance by high resolution spectroscopic imaging



# GREAT (with upGREAT & 4GREAT) has Enabled High Resolution Spectroscopic for SOFIA + Significant Imaging

 upGREAT LFA has 14 pixels at 1.9 THz and 7 pixels at 4.7 THz with limited tunability






[CII] traces component of ISM not seen in CO

The large structure to South of  $\theta^1$  Ori C is an expanding bubble. Requires high resolution imaging to understand the structure: **2600 M**<sub>sun</sub> **shell expanding at 13 kms**<sup>-1</sup>**. E** = **4x10**<sup>48</sup> **erg** 


## [CII] Line Profiles are Complex due to Source Structure and Line of Sight Absorption




Gerin+ 2015 Herschel HIFI & PACS

The low velocity resolution PACS spectra give completely erroneous picture of the [CII] emission from the source due to blending with foreground absorption by diffuse LOS clouds

Herschel HIFI




Herschel PACS R~ 1200 (250 kms<sup>-1</sup>)




What we Learn from Fine Structure Line Emission from Nearby Galaxies is Greatly Enhanced by Resolving Spectral Lines (R > 10<sup>5</sup>)

[CII] Emission from M51 (Pineda+ 2020)





# [OI] Fine Structure Line Emission has been Difficult to Interpret



Oxygen is the  $3^{rd}$  most abundant element IP of  $O^0 = 13.614$  eV, just greater than  $H^0$ 


In Photon Dominated Regions (PDRs) associated with massive young stars, oxygen fine structure line emission is extremely strong, and the 63  $\mu$ m line is widely used as a tracer of star formation

#### **QUESTIONS:**

How much neutral atomic oxygen is there in GMCs?

How does [OI] 63 µm trace star formation?







Herschel/PACS data (Gerin+ 2015)

## [OI] Mapping of W3 Galactic GMC/PDR

Goldsmith, Langer, Seo 2020





### [OI] 63 µm Self-Absorption Evident in ~1/2 of Clouds Surveyed


W3 is likely the most extreme case with absorbing  $O^0$  having  $\tau = 7$ 

Orion shows no clear evidence of selfabsorption


The radiation from O<sup>0</sup> in "hot" portion of PDR largely absorbed by intervening low-excitation material

Result is **vastly weakened** [OI] 63  $\mu$ m line Unresolved spectroscopy of 145  $\mu$ m line suggests it is relatively unaffected


What we see is highly dependent on **geometry** 




# Herschel/PACS Gave Tantalizing Glimpse of Fine Structure Line Emission from (U)LIRGS



### Tracing the Ionized ISM – Direct Probe of Massive **Star Formation**



#### NII 205 µm SOFIA/GREAT



Most are well-behaved but some have hints of complex structure that deserves follow up Why is K39 10x stronger than all the others?

### Need: High Spectral Resolution Imaging Receiver for SOFIA Studies of FIR Fine Structure & Molecular Lines

#### **FREQUENCY COVERAGE**

• [NII] 205 μm 146	1 GHz
--------------------	-------

• Multiple high-J CO lines; other FS lines

RESOLUTION R =  $10^6$  (0.3 km s<sup>-1</sup>)

#### PIXEL COUNT AND MAPPING SPEED

- Time on SOFIA is limited and expensive!
- Combine
  - State of the art noise temperature
  - High coupling efficiency
  - Extended bandwidth
  - Maximum pixel count

#### WHAT IS FEASIBLE?

1400 - 2100 GHz now!

{3300 – 5800 GHz with development}


4 GHz IF bandwidth (650 kms<sup>-1</sup> for [CII])

128 pixel square array

8X8 array times 2 polarizations


### Enabling Technology for Next-Generation Submm Heterodyne Arrays I. — REPRODUCIBLE, HIGH SENSITIVITY HEB MIXERS

**Heterodyne:** Combine (weak) astronomical signal with (strong) Local Oscillator (LO) signal in nonlinear device. Then amplify the *difference* frequency (called IF for intermediate frequency) signal generated in low noise amplifier, and finally analyze spectrally



Hot Electron Bolometer (HEB) mixers are widely used throughout submillimeter. Used in Herschel/HIFI, SOFIA/GREAT, STO2, GUSTO, ASTHROS

Well-developed Undemanding (T > 6 K) thermally Low LO power (few  $\mu$ W) 4GHz IF bandwidth  $\Leftrightarrow$  650 kms<sup>-1</sup> @ [CII]





### **HEB Array Technology**

- NbN devices routinely processed in MDL at JPL
- Individual devices diced from wafer
- Dropped into backshort/holder coupling to waveguide input
- Held in place by frontpiece consisting of 16 drilled multi-flare-angle feedhorns\*
- 4 frontpieces + backshort make 64pixel mixer array


\* These feedhorns offer broadband, high efficiency coupling to symmetric Gaussian beam

## HEB Mixer Performance is Well-Established: T(DSB) < 1000 K between 1.4 THz and 5 THz





### Array Concept: Replicate Single-Pixel Chain




Individual final LO frequency multiplication stages for EACH mixer allows OPTIMIZATION of power and achieving lowest noise temperature

Not shown: input polarizing grid separates linear polarizations. 64 pixels in 8x8 array per polarization

Very low risk; simple assembly, alignment, and cooling


# Heterodyne Focal Plane Array Imager on SOFIA - simple and low risk




Heterodyne Focal Plane Array Optics – Low

Gaussian beam telescope

Loss and Very Straightforward





### **Enabling Technology for Next-Generation Submm Heterodyne Arrays II. – WIDELY TUNABLE LOCAL OSCILLATOR SOURCE**



J. Siles, JPL

#### Huge LO Improvements since Herschel/HIFI





16-pixel LO module (~2016)

**Insight:** If the LO system produces much more power than required, it is straightforward to

- Optimize LO power for each pixel, pumping even "demanding" pixels
- Pump N > 100 pixels in array with high yield and minimal rework/replacement

100 required

Single output horn antenna

1350 1390 1430 1470 1510 1550 1590 1630 1670 1710 1750


Output Frequency (GHz)

Jet Propulsion Laboratory
California Institute of Technology

#### **Ultra-Broadband LOs**



1350 GHz -2100 GHz range



#### Ultra-Broadband LOs: First Integrated 1350-2070 GHz LO System

#### 4 pixel frequency multiplied source with two bands integrated into one single channel

One order of magnitude reduction on size, mass and power



### Enabling Technology for Next-Generation Submm Heterodyne Arrays III – DIGITAL ASIC SPECTROMETERS

**FPGA spectrometers** have worked well for GREAT. **However** – while power is modest (~20W), for large-format arrays it becomes excessive (2.5 kW)


**The solution** is system on chip digital ASIC spectrometer developed by JPL & UCLA

CURRENT: SVII 65 nm 3 GHz BW, 4096 channels 0.6 W

NEARLY DONE: SVIII – 6.5 GHz BW, 8192 channels, 4


bits, 1.7 W


Have flown on airplanes & balloons (ReckTangLE); being baselined for ASTHROS APRA balloon (launch 2023) and proposed for Astrophysics Pioneer Mission



Zhang+ 2018 IEEE Int'l. Symp. Circuits & Systems

#### **SVIII ASIC Digital Spectrometer Chip & Board**





0 - 6 GHz (or 6 - 12 GHz aliased)

**Analog inputs (ADC included)** 

Program for space qualification (radiation, temperature...) has been started, but suitable for SOFIA at the present time

## Current Status and What Should be Done for the Future of SOFIA

- Velocity resolved spectroscopy has been huge strength for SOFIA thanks to GREAT instrument. It can be even more productive in the future
- To have dramatically greater scientific impact we need to increase mapping speed by an order of magnitude and be able to sample all major components of ISM – (excited) molecular, PDR, and ionized gas
- A 128-pixel spectroscopic array will have an order of magnitude increase in mapping speed compared to upGREAT LFA and hugely more compared to HFA and 4 GREAT, and can cover needed frequency range
- This can be implemented now, but technology development program will® be valuable to
  - Improve **tuning range** of LO sources
  - Support development of frequency-multiplied sources above 2 THz
  - Increase IF **bandwidth** of mixers above 2 THz support development of MgB<sub>2</sub> mixers